

HAMAMATSU

DATA SHEET

BT(Back-thinned)-CCD Cooled Digital Camera ORCA II -BT-512




▲ Hermetic vacuum sealed air-cooled head type

The ORCA II -BT-512 features an original and unique Hamamatsu S7170 CCD chip packaged in a proprietary permanently sealed vacuum chamber evacuated to 10^{-7} Torr. This high resolution, back thinned, back illuminated CCD offers very high quantum efficiency over the spectrum from 200 nm to 1000 nm. With a huge full well capacity, low read noise and MPP (Multi-Pinned Phase) technology in the drive circuits to reduce dark current, this camera will produce rapid exposures and very high dynamic range images. Dual mode digitization offers a software selectable choice of speed or very low noise readout methods with 12 bit to 16 bit precision. Special analog contrast enhancement circuits increase versatility for even the most difficult imaging conditions.


APPLICATIONS

- Luminescence and Fluorescence
- Video microscopy
- Semiconductor imaging
- X-ray applications
- Neutron radiography
- Scintillator readout

SYSTEM CONFIGURATION



SPECTRAL RESPONSE CHARACTERISTIC



* This is typical, not guaranteed.

FEATURES

- High resolution format (512 × 512 pixels)
- High quantum efficiency from UV to NIR
- Very large full well capacity (230000 electrons typ.)
- Low readout noise design (8 electrons r.m.s. typ.)
- Software selectable dual digitizers
- Analog contrast enhancement

TYPE NUMBER

C4742-98-26L



Cooling method
A: Air-cooling
W: Water-cooling

SPECIFICATIONS


Type number	C4742-98-26LA	C4742-98-26LW
Camera head type	Hermetic vacuum sealed air-cooled head	Hermetic vacuum sealed water-cooled head
Circulating water cooler (sold separately)	-	Required
Mechanical shutter	Built-in (Control: OPEN / CLOSE / AUTO)	
Imaging device	S7170 full-frame transfer CCD	
Effective no. of pixels	512 (H) × 512 (V)	
Cell size	24 μm (H) × 24 μm (V)	
Effective area	12.29 mm (H) × 12.29 mm (V)	
Pixel clock rate	156 kHz/pixel (High-precision readout) / 2.5 MHz/pixel (High-speed readout)	
Frame rate	0.5 frame/sec (High-precision readout) / 7.0 frame/sec (High-speed readout)	
Readout noise (r.m.s.) (High-precision readout)	8 electrons (typ.) / 6 electrons (min.)	
Full well capacity	230000 electrons	
Dynamic range* (High-precision readout)	28750 : 1 (typ.) / 38333 : 1 (max.)	
Cooling method	Peltier cooling / forced-air-cooling + hermetic sealing	Peltier cooling / water-cooling + hermetic sealing
Cooling temperature	- 55 °C	- 60 °C
Dark current	0.8 electrons/pixel/s	0.3 electrons/pixel/s
A/D converter (High-precision readout)	16 bit	
Output signal (High-precision readout)	RS-422A 16 bit parallel output	
Exposure time	20 ms or more	
External control	RS-232C	
Sub-array	Yes	
External trigger	Yes	
Contrast enhancement	1, 4, 16 times (High-precision readout) / 1 to 6 times (High-speed readout)	
Lens mount	C-mount	
Line voltage	AC 100 V / AC 117 V / AC 220 V / AC 240 V, 50 Hz/60 Hz	
Power consumption	approx. 220 V·A	
Ambient storage temperature	-10 °C to + 50 °C	
Ambient operating temperature	0 °C to + 40 °C	
Ambient operating/storage humidity	70 % max. (with no condensation)	

Binning		1 × 1	2 × 2	4 × 4	8 × 8	Subarray		256 × 256	128 × 128
Frame rate	High-precision readout	6.34 frame/s	9.84 frame/s	13.6 frame/s	16.8 frame/s	Frame rate	High-precision readout	9.84 frame/s	13.6 frame/s
	High-speed readout	0.55 frame/s	1.07 frame/s	2.05 frame/s	3.75 frame/s		High-speed readout	1.86 frame/s	5.17 frame/s


* Calculated from the ratio of the full well capacity and the readout noise.

DIMENSIONAL OUTLINES (Unit: mm)


- Hermetic vacuum sealed air-cooled head (approx. 2.5 kg)




- Camera controller (approx. 8.5 kg)



- Hermetic vacuum sealed water-cooled head (approx. 2.5 kg)



- Camera controller (approx. 8.5 kg)



★ Product and software package names noted in this documentation are trademarks or registered trademarks of their respective manufacturers.

- Subject to local technical requirements and regulations, availability of products included in this promotional material may vary. Please consult with our sales office.
- Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions.
- Specifications and external appearance are subject to change without notice.

© 2001 Hamamatsu Photonics K.K.

HAMAMATSU

Homepage Address <http://www.hamamatsu.com>

HAMAMATSU PHOTONICS K.K., Systems Division
812 Joko-cho, Hamamatsu City, 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-435-1574, E-mail: export@sys.hpk.co.jp

U.S.A. and Canada: Hamamatsu Photonic Systems: 360 Foothill Road, Bridgewater, N.J. 08807-0910, U.S.A., Telephone: (1)908-231-1116, Fax: (1)908-231-0852, E-mail: usa@hamamatsu.com

Germany: Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-2658, E-mail: info@hamamatsu.de

France: Hamamatsu Photonics France S.A.R.L.: 8, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10, E-mail: infos@hamamatsu.fr

United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire, AL7 1BW, U.K., Telephone: (44) 1707-294888, Fax: (44) 1707-325777, E-mail: info@hamamatsu.co.uk

North Europe: Hamamatsu Photonics Norden AB: Smidesvägen 12, SE-171-41 Solna, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01, E-mail: info@hamamatsu.se

Italy: Hamamatsu Photonics Italia S.R.L.: Strada della Mois, 1/E 20020 Arese (Milano), Italy, Telephone: (39)02-935 81 733, Fax: (39)02-935 81 741, E-mail: info@hamamatsu.it

Cat. No. SICS1091E04
SEP/2004 HPK
Created in Japan (PDF)